



About Us:

Katena is the largest independently owned American company dedicated exclusively to ophthalmic surgical instrumentation. Katena markets its products directly to customers in the United States and through a global network of distributors in more than 110 countries around the world. Established 1975, Katena is committed to:

- Developing innovative instruments to support new surgical techniques and technological advances in the field of ophthalmology.
- Producing instruments of consistently excellent quality.
- Providing unexcelled service to its customers and strong support for its distributors worldwide.

お知らせ

このカタログの商品は、できる限り番号順に記載されています。目次は、カタログ番号の最初につく番号 (例えばK1, K2 ・・・) によって、製品のグループを特定しています。

このカタログの写真は特に注記したものを除き実物大です。イラストは細部を示す為の拡大図です。

このカタログに記載されている手術用鋼製小物は、特に注記してあるもの以外はステンレススチール製です。

すべての仕様は、予告なしに変更される場合があります。

このカタログに使用されている写真、イラストおよびテキストは著作権で保護されており、 カティーナ社の書面による同意なく複製できません。

Katena@, K-Sponge@, Performance™, Phaco-4™, Ultrakleen™, Designed for Sight™は、Katena Products社の登録商標です。

IOPatch, BioElevation, BioDome, Ambio5, keraSys and tarSys are registered trademarks of Innovative Ophthalmic Products, Inc. Ambio2 and AmbioDisk are trademarks of Innovative Ophthalmic Products, Inc. SutureGroove is a registered trademark of Advanced Meditech International – and is used with permission. Tutoplast is a registered trademark of Tutogen Medical – and is used with permission.

STAAR@, ICL™はカリフォルニア州モンロビアのSTAAR Surgical社の登録商標です。

IntraLase@はカリフォルニア州アーバインのIntraLase@社の登録商標です。



カストロビエホ氏強角膜剪刀 小刃

Castroviejo

Corneoscleral Scissors small blades **K4-1000** right (shown)

K4-1010 left



カストロビエホ氏強角膜剪刀

中刃

Castroviejo

Corneoscleral Scissors medium blades K4-1020 right (shown) K4-1030 left



トラウトマン・カストロビエホ氏強角膜剪刀

小刃

Troutman-Castroviejo

Corneoscleral Scissors with stop, small blades K4-1200 right (shown) K4-1210 left





トラウトマン・カストロビエホ氏強角膜剪刀

中刃

Troutman-Castroviejo

Corneoscleral Scissors with stop, medium blades K4-1220 right

K4-1230 left (shown)



トラウトマン・カッチン氏角膜移植剪刀

Troutman-Katzin

Corneal Transplant Scissors with stop strong curve, small blades K4-2000 right K4-2010 left (shown)



リッチ氏マイクロ角膜移植剪刀

Rich

Corneal Transplant Scissors extremely delicate, strongly curved blades For trimming the recipient corneal bed and donor button.
Can also be used for open sky anterior capsulotomy.

K4-2014 right **K4-2015** left (shown)





テリー氏マイクロ角膜剪刀

Terry

Micro Corneal Scissors round handle for meticulous control K4-2016 right (shown) K4-2017 left



ホランド氏DALK剪刀

Holland

DALK Scissors spatulated tip K4-2018 right (shown) K4-2019 left



カッチン氏角膜移植剪刀

Katzin

Corneal Transplant Scissors strong curve, medium blades K4-2020 right K4-2030 left (shown)





カストロビエホ氏角膜剪刀

Castroviejo

Universal Corneal Scissors curved, blunt tips

K4-2200 small blades **K4-2220** medium blades

K4-2240 large bladesK4-2260 strongly curved

large blades







K4-2260



カストロビエホ氏角膜移植剪刀

Castroviejo

Keratoplasty Scissors angled to side blunt tips, medium blades K4-2850



マクファーソン・ウェストコット氏結膜剪刀

McPherson-Westcott

Conjunctival Scissors curved blunt tips, small blades K4-3000





サーキシアン・ウェストコット氏剪刀

Sarkisian-Westcott

Scissors small blades with nine 1mm landmarks laser etched on one blade For intra-operative cutting and measuring with a single instrument. K4-3002



ウェストコット氏切腱剪刀

中刃

Westcott

Tenotomy Scissors blunt tips, medium blades curved K4-3004



ウェストコット氏切腱剪刀

スタンダード刃

Westcott

Tenotomy Scissors blunt tips, standard blades curved

K4-3100 right (shown)

K4-3110 left





チボー氏涙点剪刀

Tibolt

Punctal Scissors one conically shaped blade with micro grooves and textured surface For simple "Two-snip Punctoplasty". K4-3130



ウェストコット氏剪刀 ミニ

Mini Westcott

Tenotomy Scissors 4mm long curved blades blunt tips K4-3200



マクファーソン・ウェストコット氏截除剪刀

McPherson-Westcott

Stitch Scissors very sharp pointed tips small blades curved K4-4000





ウェストコット氏截除剪刀

Westcott

Stitch Scissors very sharp pointed tips standard blades curved K4-4100



ジャフィー氏截除剪刀

Jaffe

Stitch Scissors very sharp pointed tips medium blades curved K4-4150



カティーナマイクロ結紮刃

縫合糸結紮・カッティング用

Katena

Micro Tying Scissors For tying and cutting 9-0 to 11-0 sutures. 2.75mm long tying platform Not to be used as a needle driver.

K4-4200









虹彩剪刀

Iris

Scissors sharp tips, small blades K4-4400 straight

K4-4410 curved



ファイン氏截除剪刀

Fine

Stitch Scissors sharp tips, medium blades K4-4500 straight

K4-4510 curved



バラッケ氏虹彩剪刀

Barraquer

Iris Scissors blunt tips, 7mm long blades K4-4900





Capsulotomy Scissors

切囊剪刀

Capsulotomy Scissors

バンナス氏切嚢剪刀

Vannas

5mm long blades, sharp tips

K4-5000 straight

K4-5010 curved

K4-5020 angled forward

K4-5030 angled to side

ギルス・バンナス氏切嚢剪刀

Gills-Vannas

7mm long blades, sharp tips

K4-5050 straight

K4-5060 curved

K4-5070 angled forward

ホランド氏DALK剪刀

Holland DALK Scissors

4.5mm angled blades

K4-5075 blunted tips

マクファーソン・パンナス氏虹彩剪刀

McPherson-Vannas

8mm long blades, sharp tips

K4-5100 curved











K4-5010

K4-5020

K4-5030

K4-5050









K4-5075 K4-5100

切囊剪刀

Capsulotomy Scissors

Katena-Vannas

7mm extra thin long blades, sharp tips

カティーナ・パンナス氏切嚢剪刀

K4-5110 Katena-Vannas, straight

K4-5111 Katena-Vannas, curved

ギルス・バンナス氏切嚢剪刀

K4-5112 Gills-Vannas, angled forward 10mm extra thin long blades, sharp tips

カティーナ・パンナス氏切嚢剪刀

K4-5120 Katena-Vannas, straight

K4-5121 Katena-Vannas, curved

スターン・ギルス氏切嚢剪刀

K4-5122 Stern-Gills, angled forward

ウリベ・スターン氏切嚢剪刀

K4-5123 Uribe-Stern, curved,

angled, forward

10mm long angled blades, curved tips

ギルス氏切嚢剪刀

K4-5125 Gills, sharp points

ウトホフ・ギルス氏切嚢剪刀 K4-5126 Uthoff-Gills, blunt points











K4-5111

K4-5112

K4-5120

K4-5121





K4-5122

K4-5123

K4-5125

K4-5126



硝子体剪刀

20 G

Vitrectomy Scissors

20 gauge shank Stainless steel construction with serrated plastic handle. Autoclavable.



Scissors

horizontal opening straight K4-5200 ホリゾンタル 直



Scissors vertical opening angled

60° K4-5240 バーチカル



Scissors

horizontal opening angled 30° K4-5210 ホリゾンタル # 30°



Scissors vertical opening angled

K4-5250 バーチカル ## 80°

80°



硝子体・網膜用剪刀 チタニウム製スクウィーズハンドル Katena

"Squeeze Handle" Scissors These "squeeze handle" instruments feature titanium handles to reduce overall weight, and front ends of hardened stainless steel for secure cutting.



曲 60°

Scissors

horizontal, straight **K4-5300** 20ga. ホリゾンタル



Scissors

vertical, 40° K4-5330 20ga. バーチカル、40°



Scissors

horizontal, curved **K4-5310** 20ga. K4-5423 23ga. K4-5425 25ga. ホリゾンタル、反



Scissors

vertical **K4-5433** 23ga. バーチカル



Scissors

horizontal, angled 45° **K4-5320** 20ga. ホリゾンタル 曲 45°



Scissors

vertical, 80° K4-5340 20ga. バーチカル、80°



アリオ氏MICS剪刀

Alio

MICS Scissors
23 gauge curved shaft with
horizontal micro blades
For cutting delicate membranes.

K4-5351





シメカ氏IOLループ切断剪刀

Scimeca

Lens Loop Amputator For cutting polypropylene or PMMA lens haptics in the anterior or posterior chamber.

K4-5550 with guard



IOL Cutters

ドネンフェルド氏マイクロIOLカッター

Donnenfeld

Micro IOL Cutter 4mm long blades micro notches in blades for traction For cutting soft IOLs through a 1.5mm incision. K4-5560




Used with K5-7570

ストルテ氏ソフトIOLカッター チタニウム製

Stolte

IOL Cutter titanium slotted cover to retain soft IOL during cutting process For cutting soft IOLs.

K4-5565



ソフトIOLカッター

Katena

IOL Cutter 6mm long serrated blades For cutting soft IOLs.

K4-5571

