



katena[®] OPTHALMIC
INSTRUMENTS

Reusable Ophthalmic Instruments

K5 鑷子等

OPHTHALMOLOGY

corzaeye.com

 **corzamedical**

Committed to developing innovative instruments that support new surgical techniques and technological advancements in ophthalmology.

Katena markets its products directly to customers in the United States and through a global network of distributors in over 110 countries.

Table of Contents

Catalog Number Prefix	Product Group	Page Number
K1	Speculums • Scleral Depressors • Lid Plates • Fascia Needle • Retractors	4 - 21
K2	Knives • Spuds • Burrs • Blade Breakers • Diamond Knives Trabeculectomy Punches	22 - 30
K3	Spatulas • Choppers • Manipulators • LASIK Spatulas • Lens Loops Manual Nucleus Fragmentation • Iris Retractors • Hooks • Fixation Rings Muscle Hooks • Corneal Markers • Toric IOL Markers • Calipers	31 - 80
K4	Scissors • IOL Cutters	81 - 95
K5	Forceps	96 - 130
K6	Needle Holders	131 - 136
K7	Rongeurs • Lacrimal Instruments • Cannulas Irrigating Manipulators & Choppers • I/A Systems	137 - 165
K8	Eye Spheres • Bipolar Forceps	166 - 167
K9	Sterilizing Cases • Eye Shields • Accessories	168 - 175
K20	Disposables - Trephine Blades • Blades • Sponges	176 - 179
	Index	180 - 183

Forceps

Gills-Colibri

Tissue Forceps
very fine pointed tips
with tying platform

K5-I450



Colibri Style

Tissue Forceps
very delicate, 1x2 teeth, 0.12mm

K5-I500 Colibri

K5-I520 Harms-Colibri, with tying platform

KT5-I520 Harms-Colibri, titanium, with tying platform

K5-I522 Dyson-Colibri, short tips, extremely thin shanks, with tying platform

Standard, 1x2 teeth with tying platform

K5-I550 Colibri, 0.4mm

Corneal Suturing Forceps

very delicate, 1mm spread, 1x2 teeth, 0.12mm

K5-I560 Polack

Corneal Forceps

very delicate, 1x1 teeth, 0.1mm

K5-I700 Pierse-Colibri



K5-I500



K5-I520




K5-I522



K5-I550



K5-I560



K5-I700

Forceps

Hofmann-Polack

Corneal Suturing Forceps
with guide pin
1mm spread
1x1 teeth, 0.2mm


K5-1566



Hoskin

Tissue Forceps, curved

K5-1716 #16, fine grooved tips



Hoskin

Tissue Forceps, straight


K5-1718 #18, fine grooved tips

K5-1719 #19, micro grooved tips


K5-1728 #28, very delicate,
micro grooved tips



K5-1718



K5-1719



K5-1728



Forceps

Hoskin


Tissue Forceps, Colibri style

K5-1722 #22, micro grooved tips

K5-1730 #30, very delicate
micro grooved tips



K5-1722



K5-1730



Pierse

Corneal Forceps

very delicate

1x1 teeth, 0.1mm

straight

K5-1800



Fechtner


Conjunctiva Forceps

very delicate

ring shaped tips

with tying platform

K5-1820



Castroviejo-Colibri

Forceps

very delicate

1x2 teeth, 0.12mm

K5-1900



Forceps

Harms-Colibri

Forceps
very delicate
1x2 teeth, 0.12mm
with tying platform

K5-1902



Hyde

Corneal Forceps
2x2 teeth, 0.12mm
with tying platform

K5-1920



Bonn

Forceps
very delicate
1x2 teeth, 0.12mm

K5-2000



Bonn

Forceps
1x2 teeth, 0.12mm
with tying platform

K5-2020 stainless steel

KT5-2020 titanium




Forceps

Thorpe

Corneal Forceps
2x3 teeth, 0.12mm
with tying platform

K5-2100




McPherson

Tissue Forceps
1x2 teeth, 0.4mm
with tying platform

K5-2300 straight

K5-2310 angled



Kelman-McPherson

Forceps
8mm long jaws
1x2 teeth, 0.3mm

K5-2330



Forceps

Castroviejo

Suturing Forceps
1x2 teeth
with tying platform

K5-2500 0.12mm

K5-2510 0.3mm

K5-2520 0.5mm

K5-2530 0.9mm

KT5-2500 0.12mm, titanium



Snider-Castroviejo

Suturing Forceps
with guide pin
1x2 teeth, 0.12mm
with tying platform

K5-2550



Stern-Castroviejo

Fixation Forceps
1x2 teeth, 0.5mm
with lock

K5-2552




Moody

Fixation Forceps
1x2 teeth, 0.5mm
with lock

K5-2553 left (not shown)

K5-2554 right




Forceps

Helveston

Fixation Forceps
1x2 teeth, 0.6mm

*For globe fixation
during strabismus
procedures.*

K5-2556



Girard

Corneoscleral Forceps
1x2 teeth, 0.12mm
with tying platform

K5-2900 for use in right hand

K5-2910 for use in left hand



Perl

Corneal Forceps
1x2, 0.08mm teeth

*For atraumatic grasping
and fixation of tissue.*

K5-2914 for use in left hand



Maumenee-Colibri

Corneal Forceps
very delicate
1x2 teeth, 0.12mm

K5-3000




Forceps

Bishop-Harmon

Tissue Forceps
1x2 teeth

K5-3500 delicate


K5-3510 standard (not shown)



Bishop-Harmon

Dressing Forceps
criss-cross serrated tips

K5-3550



Fixation Forceps

1x2 teeth, 0.3mm

K5-3560 Sauer

K5-3562 Manhattan Eye & Ear,
with tying platform



Dressing Forceps

Serrated, delicate

K5-4000 straight

K5-4010 curved



Forceps

Tissue Forceps

1x2 teeth

K5-4100 delicate, straight

K5-4110 delicate, curved (not shown)

K5-4140 standard, straight (not shown)



Adson


straight, 5" (125mm) long

K5-4150 Dressing Forceps
serrated jaws,
delicate

K5-4152 Dressing Forceps
serrated jaws,
standard

K5-4170 Tissue Forceps
1x2 teeth, delicate

K5-4172 Tissue Forceps,
1x2 teeth, standard



1x2 teeth, delicate



Bayonet

Dressing Forceps
serrated jaws
5 3/4" (146mm) long

K5-4180



Lester

Fixation Forceps
straight

K5-4300 1x2 teeth

K5-4320 2x3 teeth (not shown)



Forceps

Fixation Forceps

1x2 teeth

K5-4500 Bracken, delicate teeth

K5-4550 Elschmig, standard teeth
(not shown)



Troutman

Superior Rectus
Forceps
1x2 teeth
angled

K5-4700



Tying Forceps

For 8-0 to 11-0 sutures.

McPherson

5mm long smooth jaws

K5-5000 straight

KT5-5000 straight, titanium

K5-5010 angled

KT5-5010 angled, titanium

Kelman-McPherson

7.5mm long smooth jaws

K5-5030 angled

KT5-5030 angled, titanium

Bechert


10mm long smooth jaws

K5-5049 angled


Gaskin

11mm long smooth jaws


K5-5050 angled




K5-5000




K5-5010



K5-5030



K5-5049



K5-5050

Forceps

Haberle-McPherson


Lens Holding Forceps

7.5mm long jaws

with groove, angled

For grasping IOL haptic.

K5-5040



Kurstin

Flap Stretching Forceps


disc shaped

serrated jaws

Designed for stretching

the flap to eliminate striae.

K5-5056



Buratto

LASIK Flap Forceps

disc shaped

serrated jaws


K5-5058



Mendez SMILE Forceps

For removing the lenticule.

K5-5062



Forceps


Inamura

Capsulorrhexis Forceps
delicate tapered shanks
to seal the incision during
continuous curvilinear
capsulorrhexis

Excellent for MICS.

K5-5070 short tips

K5-5072 long tips



Utrata

Capsulorrhexis Forceps
very delicate, short handle
triangular grasping tips
extremely thin with landmarks
at 2.5mm and 5.0mm

K5-5071 straight



Masket

Capsulorrhexis Forceps
very delicate, short handle
triangular grasping tips
extremely thin with landmarks
at 2.5mm and 5.0mm

K5-5074 vaulted



Utrata

Capsulorrhexis Forceps Short
Handle, very delicate, triangular
grasping tips

K5-5075 3-hole handle


K5-5076 serrated handle



Kraff

Capsule Tag Forceps
11mm long shanks
1.5mm long grasping tips

K5-5080



Forceps

Utrata

Capsulorrhexis Forceps
very delicate
triangular grasping tips
extremely thin
11mm long shanks

K5-5081 stainless steel

KT5-5081 titanium

Also available for MICS, see K5-5091.



Utrata

Capsulorrhexis Forceps
very delicate triangular
grasping tips extremely
thin 11mm long shanks
round handle with guide
pin for alignment

K5-5082




Nevyas

Capsulorrhexis Forceps
cystotome shaped tips
extremely thin
11mm long shanks

*For tearing and grasping
the anterior capsule.*

K5-5083

Also available for MICS, see K5-5093.



Masket

Capsulorrhexis Forceps
very delicate
triangular grasping tips
extremely thin 11mm long
vaulted shanks

K5-5084



Forceps

Masket

Capsulorrhexis Forceps
very delicate
triangular grasping tips
extremely thin 11mm long
vaulted shanks
round handle with guide pin
for alignment

K5-5086



Jaffe

Capsulorrhexis Forceps
very delicate
round grasping tips
extremely thin
11mm long shanks

K5-5088




Giannetti

MICS Capsulorrhexis Forceps
ultra delicate shanks and tips
maximum jaw opening 1.75mm

*For performing capsulorrhexis
through a micro incision.*

K5-5090



Utrata

MICS Capsulorrhexis Forceps
very delicate
triangular grasping tips
thin 11mm shanks
maximum jaw opening 1.75mm

*For performing capsulorrhexis
through a micro incision.*

K5-5091



Forceps

Giannetti

MICS Capsulorrhexis Forceps
ultra-delicate shanks and tips
maximum jaw opening of
1.75mm flat handle

*For performing capsulorrhexis
through a micro incision.*

K5-5092



Nevyas

MICS Capsulorrhexis Forceps
cystotome shaped tips
ultra-delicate shanks and tips
maximum jaw opening of 1.75mm
flat handle

*For penetrating and grasping the anterior
capsule while performing a Capsulorrhexis
through a micro incision.*

K5-5093



Florakis

Endothelial Forceps
extra-delicate, reversed
triangular tips angled 75°

*To aid in the removal of endothelium
from the host cornea.*

K5-5095



Anis

Tying Forceps
extra delicate smooth jaws

For 9-0 to 11-0 sutures.

K5-5113 curved



Forceps


Jaffe

Tying Forceps
very delicate smooth jaws
6mm long platform

For 8-0 to 11-0 sutures.

K5-5200 straight

K5-5210 curved



Gerl

Tying Forceps
curved shanks, 6mm long platform

For 9-0 to 11-0 sutures.

KT5-5215



Tennant

Tying Forceps
extra delicate smooth jaws with 6mm long
platform round handle with guide pin


For 9-0 to 11-0 sutures.

K5-5220 straight

KT5-5220 straight, titanium

K5-5230 curved

KT5-5230 curved, titanium



Harms

Tying Forceps
standard, smooth jaws

For 7-0 to 9-0 sutures.

K5-5400 straight

K5-5410 curved



Helveston

Tying Forceps, standard, smooth jaws

For 4-0 to 6-0 sutures.

K5-5420



Forceps

Fine

Tying Forceps
standard, smooth jaws
For 7-0 to 9-0 sutures.

K5-5550



Arita

Meibomian Gland Compressor
broad, smoothly polished jaws

K5-5900



Barraquer

Cilia Forceps
lightly angled
Textured jaws provide superior grasping capabilities.

K5-6000



Cilia Forceps

slanted jaws
gold plated

K5-6020



Gradle

Cilia Forceps
smooth jaws

K5-6100




Forceps

Bonaccolto

Conjunctiva Forceps
longitudinally serrated jaws
with cross serrations at tips

K5-6200 1.2mm wide

K5-6210 1.7mm wide (not shown)



Nugent

Utility Forceps

K5-6214 smooth jaws (not shown)

K5-6215 serrated jaws



Gold

Punctal Plug Forceps
0.2mm wide longitudinal
groove inside jaws
pointed tips

K5-6218



Hersh-Greenstein-Gelles

CTAK Inlay Forceps
Curved, smooth tissue forceps for
placement of the CTAK inlay into
the fs laser channel.

K5-6220



Forceps

Foreign Body

Forceps
straight

K5-6400



Utility & Splinter

Forceps
pointed, short jaws
Suitable for epilation.

K5-6500 #H



Jewelers

Forceps
standard, straight

K5-6510 #1



Jewelers

Forceps
delicate, straight

K5-6520 #3-C



Jewelers

Forceps
delicate, straight

K5-6530 #3



Forceps

Jewelers

Forceps
delicate, straight

K5-6540 #4



Jewelers

Forceps
extra delicate, straight

K5-6550 #5



Jewelers

Forceps
extra delicate, angled

K5-6560 #5/45



Jewelers

Forceps
curved

K5-6570 #7



Suture

Removal Forceps
very fine pointed tips
with textured jaws

Ideal for grasping wet sutures.

K5-6575



Forceps

Arruga

Capsule Forceps

K5-7000



Akahoshi

Phaco PreChopper
4.5mm long delicate jaws
with sharp tips

*For cracking, dividing and rotating
the nucleus without sculpting.*

K5-7230




Akahoshi

Combo PreChopper

*For cracking, dividing and rotating
both hard and soft nuclei without sculpting.*

K5-7232



Forceps

Ernest

Nucleus Cracker
extra delicate cross-action
serrated paddle tips

K5-7240 For cracking the nucleus
through a 2.5mm incision.

K5-7242 For cracking the nucleus
through a 2mm incision.




Dodick

Nucleus Cracker
cross-action
longitudinally serrated
paddle shaped jaws

*For cracking the nucleus
through a 3mm incision.*

K5-7244



Hu

Femto Forceps
blunted 0.3mm tips
cross-action mechanism

*For opening incisions
created with a
femtosecond laser.*

K5-7245



Alfonso

Nucleus Grasping Forceps
two short rows of
delicate interlocking teeth

*For grasping and extracting
nucleus fragments.*

K5-7253




Forceps

Femto


Nucleus Forceps
22 gauge
flattened textured tips
angled at 90°
reverse action
*For removing a complete cap,
finishing a rhexis and
splitting nucleus fragments.*

K5-7270




Vitreo - Retinal Forceps


Stainless steel with
serrated plastic handle
autoclavable




smooth jaws
20 gauge, straight
K5-7500




serrated jaws
20 gauge, straight
K5-7510



smooth jaws, 1x2 teeth
20 gauge, straight
K5-7520



smooth jaws
end gripping
20 gauge, straight
K5-7530




smooth jaws
20 gauge, vertical
K5-7540

Forceps

Endothelial Grasping Forceps

curved, 23 gauge
delicate end-gripping platforms
jaws 90° off-set from
“squeeze handle”
mechanism


K5-7550



ICL™

Lens Cartridge Loading Forceps
straight, 20 gauge
jaws feature a series
of concave and
convex surfaces

K5-7560



Donnenfeld

Micro IOL Holding Forceps
21 gauge blunt, paddle-shaped tips
central hole for traction and visibility

*For grasping a soft IOL during cutting
and explantation.*

K5-7570

(use with K4-5560)



Forceps


MICS Capsulorhexis Forceps

23 gauge curved shaft with short micro-grasping tips


Designed for use through a very small paracentesis.

K5-7651 Alio

K5-7655 Fine-Ikeda



K5-7651




K5-7655



Squeeze Handle Forceps


These “squeeze handle” instruments feature titanium handles to reduce overall weight, and front ends of hardened stainless steel for secure grasping.



horizontal opening
serrated jaws

K5-7600 20ga.

K5-7703 23ga.




horizontal opening
end gripping jaws

K5-7610 20ga.


K5-7713 23ga.

K5-7715 25ga.




vertical opening
smooth side
gripping jaws

K5-7620 20ga.



membrane peeling
smooth 45°
angled jaws

K5-7630 20ga.



membrane peeling
grooved 45°
angled jaws


K5-7631 20ga.



asymmetrical

K5-7723 23ga.

K5-7725 25ga.



ILM

K5-7733 23ga.

K5-7735 25ga.


Forceps

Prasad

Micro Tying Forceps
23 gauge
straight 2.65mm long jaws
with textured platforms
25mm long curved shaft

K5-7743 straight

K5-7753 angled



Watzke

Sleeve Spreading Forceps


K5-7950



Rosenwasser

Donor Lamella Inserting Forceps
thin vaulted 12mm long shanks
precalibrated heel stop to
prevent tips from closing

K5-8010



Forceps

Lens

Loading Forceps

K5-8100



ICL™

Loading Forceps

11mm long

smoothly polished jaws

K5-8180



ICL™

Removal Forceps

11mm long

finely textured jaws

K5-8181



Dodick

IOL Folding Forceps

For soft IOLs.

K5-8210




Forceps

Livernois

IOL Folding Forceps
concave jaw surfaces

For soft IOLs.

K5-8215



Steinert

IOL Inserting Forceps
conically shaped jaws
with flat inside surfaces

*Used for folding and
inserting soft IOLs.*

K5-8219




Livernois

IOL Holding Forceps
highly polished round jaws

*For holding a soft IOL
while folding with an
inserting forceps.*

K5-8221 angled



Forceps

Ernest-McDonald

IOL Inserting Forceps
5mm long straight jaws
low profile
patented cross-action

For inserting silicone IOLs.

K5-8228



Faulkner

Lens Holding Forceps
round, smooth jaws

For holding soft IOLs.

K5-8230



Buratto

Soft IOL Inserting Forceps
extra delicate highly
polished convex jaws

*For foldable acrylic and
silicone IOLs.*

K5-8268



Spaleck

Soft IOL Inserting Forceps
delicate highly polished
convex jaws

*For foldable acrylic and
silicone IOLs.*

K5-8270



Forceps

Blaydes

Lens Holding Forceps
very delicate narrow jaws
angled 45°

K5-8500



Shepard

Lens Holding Forceps
duckbill shaped jaws
gently curved

K5-8600



Forceps

Jameson

Muscle Forceps
with slide lock

K5-9000 4 teeth, child size, right

K5-9010 4 teeth, child size, left (not shown)

K5-9020 6 teeth, adult size, right (not shown)

K5-9030 6 teeth, adult size, left



Berke

Ptosis Forceps
with slide lock
longitudinally grooved jaws

K5-9300 20mm long

K5-9310 27mm long



Forceps

Ptosis

Lid Clamp
6-pin, concave
longitudinal grooves
in tips, sliding lock

K5-9315



Erhardt

Lid Forceps

K5-9400




Downes

Lid Clamp

*Designed for protecting the
eye and holding the lid margin
during upper lid laser
blepharoplasty.*

K5-9410



Forceps

Chalazion Forceps

Round

K5-9450 Lambert, 8mm I.D.

K5-9500 Hunt, 12mm I.D.




Chalazion Forceps

Oval


K5-9550 Baird, 11 x 8mm I.D.
small

K5-9600 Francis, 14 x 11mm I.D.
medium

K5-9650 Desmarres, 20 x 13mm I.D.
large



K5-9600



K5-9650

Wies

Chalazion Forceps
grooved upper jaw

K5-9750



Forceps

Serrefine

Serrated jaws
small 1½" (38mm) long

K5-9850



Hartmann

Mosquito Forceps
extra delicate
3¾" (95mm) long

K5-9900 straight

K5-9910 curved



Halsted

Mosquito Forceps
delicate
5¾" (145mm) long

K5-9940 straight

K5-9950 curved



Forceps

Baby Jones

Towel Clamp
cross-action
2¼" (55mm) long

K5-9960



Backhaus

Towel Clamp

K5-9964 3½" (90mm) long

K5-9965 5¼" (135mm) long

